




Do not anodize o-ring grooves or any surface interior to o-ring grooves. Polish interior surfaces.



ITEM	PART NUMBER	DESCRIPTION	QTY.
1	PFS15102	CCD Head	1
<b>THE OBSERVATORIES OF THE CARNEGIE INSTITUTION FOR SCIENCE</b> 813 Santa Barbara Street Pasadena, CA 91101			
MATERIAL	Al 6061-T6	UNLESS OTHERWISE SPECIFIED FRACTIONAL DECIMAL ANGULAR ± 1/32 .XX ± 0.01 ± 0.1 DEG	
EST WEIGHT	7.477	.XXX ± 0.002 GOAL .XXX ± 0.005 REQD BREAK SHARP EDGES. MACHINED SURFACES ALL DIMENSIONS IN MM	
SCALE	1:2	<b>P</b> lanet <b>F</b> inder <b>S</b> pectrograph	CCD Head
DRWN	CRANE	4/12/2017	SHEET B DRAWING PFS15102
APVD			SHT 1 OF 2 REV



THE OBSERVATORIES OF THE CARNEGIE INSTITUTION FOR SCIENCE 813 Santa Barbara Street Pasadena, CA 91101					
MATERIAL	Al 6061-T6		UNLESS OTHERWISE SPECIFIED FRACTIONAL DECIMAL ANGULAR $\pm 1/32$ .XX $\pm$ 0.01 $\pm$ 0.1 DEG		
EST WEIGHT	7.477		.XXX $\pm$ 0.002 GOAL .XXX $\pm$ 0.005 REQD BREAK SHARP EDGES. MACHINED SURFACES ALL DIMENSIONS IN MM		
SCALE	1:2		Planet Finder Spectrograph		CCD Head
DRWN	CRANE	4/12/2017	SHEET	DRAWING	REV
APVD			B	PFS15102	SHT 2 OF 2